



TecnoSUC  absolut vacuum



ANNIVERSARY



Tecno suc air

Model	Code	Dental units	Lt/min	Vuoto m bar	kW	dB	Kg	Dimensioni AxBxC mm
tecno suc air 1 carenata	GA002052	1	1200	100	0.37	54	20	270x290x330
tecno air 1	GA002053	1	1200	100	0.37	58	10	250x250x250
tecno air 2	GA002054	2	2300	150	0.75	60	14	290x270x300
tecno air 3	GA002055	3/4	3500	180	1.5	66	21	330x280x350
tecno air 5	GA002039	5/6	4350	230	2.2	66	30	380x330x380



Tecno suc water

Model	Code	Dental units	Lt/min	Vuoto m bar	kW	dB	Kg	Dimensioni AxBxC mm
tecno suc 2 water carenata	GA002056	2	383	532	0.75	60	40	450x340x460
tecno suc 2 water	GA002057	2	383	532	0.75	62	30	320x270x370
tecno suc 4 water	GA002058	4	750	560	1.5	63	40	400x320x420
eco tecno suc 2 carenata	GA002065	2	383	40	0.75	64	60	495x530x555
eco tecno suc 4 carenata	GA002066	4	750	42	1.5	64	70	495x530x640


tecno suc eco separatore di amalgama 98%



Tecno suc wet

Model	Code	Dental units	Lt/min	Vuoto m bar	kW	dB	Kg	Dimensioni AxBxC cm
VS 1 *	GA004298	1	-	160 170	0.38	63	14	260x260x360
VS 2 *	GA004299	2/3	-	170	1.1	63	25	420x410x500
VS 3 *	GA004998	3/5	-	170	1.48	64	35	420x410x550

* carenatura optional



dental units and not only